



Part **Basket without Lock**  
 Part Number **02830**  
 Serial Number **SR N11012696 – SR N11012846**



# Test Certificate

		<b>USP / Ph.Eur. Tolerances</b>		
<b>Dimensions</b>	Basket	Basket rim O.D. Ø	25.0 mm	± 3.0 mm
		Height	37.0 mm	± 3.0 mm
		Wire Ø	0.25 mm	0.22 – 0.31mm
		Wire openings	0.4 mm	0.36 - 0.44 mm
		Screen opening height	27.0 mm	± 1.0 mm
		Screen O.D. Ø	22.2 mm	± 1.0 mm
		Clear opening Ø	20.2 mm	± 1.0 mm
<b>Material</b>	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2022-01-03

Bruno Balli  
 SOTAX Quality Management  
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland