




Test Certificate

Part 1 Liter Vessel
 Part Number 14143-01
 Serial Number 84176 - 84596



		USP / Ph.Eur. Tolerances	
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm
		Height b	200 mm
		With hemispherical bottom	160 – 210 mm
Material	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2021-01-05

Bruno Balli
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland