




Part 1 Liter Vessel Inactinic Brown  
 Part Number 14143-12  
 Serial Number P36 000068 – P36 0000134



**USP / Ph.Eur.  
Tolerances**

<b>Dimensions</b>	Vessel	InsideØ a	102.9 to 105.1 mm	98 – 106 mm
		Height b	168.2 mm	160 – 210 mm
		With hemispherical bottom		
<b>Material</b>	Borosilicate glass inactinic brown			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Navi Mumbai, 2021-12-22

12/22/2021 3:58:09 PM

Moreshwar Pulekar  
 SOTAX Quality Management  
 C-64 | First Floor | TTC Industrial area | Turbhe | Navi Mumbai 40070 | India

# Test Certificate