





Part Automated Disintegration Disk Set (3 x Disk)
 Part Number 23444-01



				Ph.Eur. Tolerances
Dimensions	A	Disk outside Ø	31.40 mm	± 0.13 mm
	B	Thickness	15.30 mm	± 0.15 mm
	C	7 Holes Ø	3.15 mm	± 0.1 mm
	D	Internal bolt circle	8.40 mm	± 0.1 mm
Specific gravity			ζ 1.18 – 1.20 kg/dm ³	--
Material		PMMA		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2021-03-10

Bruno Balli
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

Test Certificate